## Walkers by Orliman®

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## Walkers by Orliman®

### INSTRUCTIONS FOR USE, STORAGE AND WARRANTY

Dear Customer,

Thank you very much for placing your trust in an Orliman product. Please read the instructions carefully. Keep these instructions and the packaging for future reference. If you have any questions or concerns, please contact your doctor, orthopaedic specialist or our customer service department.

ORLIMAN S.L.U. guarantees all its products as long as the original configuration has not been manipulated or altered except for the intended use as described in these instructions.

If the products are used in combination with other products, replacement parts or systems, make sure they are compatible and made by Orliman®. It does not guarantee any products with altered characteristics due to improper use, defects or breakage of any kind. The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. If any serious incidents related to the product occur, notify Orliman S.L.U. and the corresponding competent authority in your country.

Orliman would like to thank you for choosing this product and hopes you a speedy recovery.

### REGULATIONS

**MD** This article is defined as a class I medical device. A Risk Analysis (UNE EN ISO 14971) has been carried out, minimising the existing risks. Tests have been in accordance with European Regulation UNE-EN ISO 22523 on Prostheses and Orthoses.

### INDICATIONS

Protection and control following traumatism or surgical intervention, sprains and ligament injuries, as a method to replace plaster in the treatment of preulcerous or ulcerous conditions of the foot's plantar surface. Acute plantar fasciitis and in those cases which require a uniform distribution of body weight on the plantar surface and immobilisation of the tibia-tarsus articulation.

### FITTING INSTRUCTIONS

For best therapeutic results considering different pathologies and to extend the useful life of the product, it is essential to choose the correct size for each patient or user. Excessive compression may be intolerable; adjusting the compression to be firm yet comfortable is recommended. If the product must be adapted, any such adjustments must be done by an orthopaedic specialist or healthcare professional legally certified to do so who must make sure the end user or person responsible for fitting the product properly understands how it works and should be used.

**When fitting the product, you must adhere to the following instructions:**

#### 1 Adapting the walker to the patient:

Placing the walker in a sitting position and putting a sock (cotton, if possible) on the patient's leg and foot is recommended.

1-Separate the rigid structure from the wrap-around sleeve. Place the wrap-around sleeve snug against the patient's foot and leg and tighten it as much as possible using the quick-close system. Position the protective toe pad. If necessary, place the protective pads in sensitive areas such as the instep or malleoli.

2-Place the leg over the brace and make sure the foot is in direct contact with the base and back and that the rigid structure is centred with the leg and ankle load line.

**For Ref.: EST-083, EST-086, EST-087 and EST-088:**

3-If necessary, shape the splints for appropriate anatomic adaptation. Remove the protective plastic from the mediolateral splints to secure them to the wrap-around sleeve. For Ref. EST-087, secure the tibia protector plate in the front.

4-Close and adjust the brace with quick-close straps in the foot area.

5-Close the quick-close straps on the leg. For Ref. EST-083 and EST-086, use the straps for the tibia and calf protector plates.

6-Check for a correct fit. The product must be firmly adjusted yet comfortable.

**For Ref.: EST-089 and EST-092:**

3-Secure the tibia protector plate to the front.

4-Close and adjust the brace with quick-close straps in the foot area.

5-Close the quick-close straps on the leg.

6-Check for a correct fit. The product must be firmly adjusted yet comfortable.

### Flexion-extension control adjustment (EST-083 and EST-086):

- 1-Slide and open the cover to expose the pins.
- 2-Release the pins, select the desired position and reinsert the pins. If necessary, first loosen the locking screw to adjust the flexion-extension.
- 3-Close the cover by sliding the protective part. To lock the walker boot at a fixed angle of 90°, position both pins at 0°. When you wish to lock it in a fixed position at 0°, 7,5°, 15° or 22.5° plantar or dorsiflexion, tighten the Allen head screws on both joints.

### Using the inflation pump (EST-083, EST-089 and EST-092):

- 1-To inflate, turn the valve towards INFLATE (clockwise). The valve must be really tight so no air is lost.
- 2-Press the orange circular part several times until you get the desired pressure.
- 3-Check for a correct fit. The product must be firmly adjusted yet comfortable.
- 4-To release the air from the bag, turn the valve towards DEFLATE (counter-clockwise) until comfortable.

### PRECAUTIONS

Before each use, check that all product components are present as per the fitting process. Periodically check the conditions of the product. If you observe any defect or anomaly, immediately report it to the issuing establishment.

This product is made of inflammable material. Do not expose the products to situations that could set them on fire. In the event of a fire, quickly get them off your body and use the proper resources to extinguish the fire.

To avoid minor discomfort caused by sweating, we recommend using some type of cotton fabric to separate the skin from contact with the product material. For discomfort such as chafing, irritation and swelling, remove the product and see a doctor or orthopaedic specialist. The product should only be used on healthy skin. It is not recommended for use over open scars with swelling, redness or hotspots.

Products marked with the  symbol contain natural rubber latex and can cause allergic reactions in people sensitive to latex.

Products marked with the  symbol contain ferromagnetic components and, therefore, extreme precaution must be taken if you undergo an MRI scan or are exposed to radiation associated with diagnostic or therapeutic procedures.


### RECOMMENDATIONS-WARNINGS

The use of these products is conditioned by the indications. Although the product is not defined as a single-use device, using it on a single patient only is recommended and only for the intended purposes as described in these instructions or by a healthcare professional.

When disposing of the product and its packaging, you must strictly adhere to the legal regulations in your community.

### RECOMMENDATIONS FOR STORAGE AND WASHING

When not using the product, store it in the original packaging in a dry place at room temperature. Stick the Velcro to each other (if the orthotic device has them), frequently wash by hand with warm water (30° C max.) and mild soap. To dry the product, use a dry towel to absorb as much moisture as possible and let it dry at room temperature. Do not hang it up or iron the product and do not expose it to direct heat sources such as stoves, dryers, direct sun exposure, etc. When using or cleaning the product, do not use abrasive or corrosive substances, alcohol, ointments or liquid solvents. If not dried off properly, the detergent residue may irritate the skin and cause the product to deteriorate.

**EN** The manufacturing date is included in the batch number which can be found on the packaging label as **(00)**, in the following way: the second and third digits represent the year of manufacture and the fourth and fifth digits represent the month.

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